

Product datasheet for **GA113739**

Human C1orf198 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C1orf198 CRISPRa kit - CRISPR gene activation of human chromosome 1 open reading frame 198
Symbol:	C1orf198
Components:	GA113739G1, C1orf198 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113739G2, C1orf198 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113739G3, C1orf198 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H425</u>
Synonyms:	DKFZp667D152; FLJ14525; FLJ16283; FLJ38847; MGC10710
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)