

Product datasheet for **GA113784**

Human PWWP3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PWWP3A CRISPRa kit - CRISPR gene activation of human PWWP domain containing 3A, DNA repair factor
Symbol:	PWWP3A
UniProt ID:	Q2TAK8
Synonyms:	EXPAND1; HSPC211; MUM-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)