

Product datasheet for **GA113787**

Human WDR73 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR73 CRISPRa kit - CRISPR gene activation of human WD repeat domain 73
Symbol:	WDR73
Components:	GA113787G1, WDR73 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113787G2, WDR73 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113787G3, WDR73 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P4I2</u>
Synonyms:	GAMOS; HSPC264
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)