

Product datasheet for **GA113791**

Human SERAC1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SERAC1 CRISPRa kit - CRISPR gene activation of human serine active site containing 1
Symbol:	SERAC1
Components:	GA113791G1, SERAC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113791G2, SERAC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113791G3, SERAC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96JX3
Synonyms:	FLJ14917; FLJ30544
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)