


## Product datasheet for **GA113802**

### Human KIAA1984 (CCDC183) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | CCDC183 CRISPRa kit - CRISPR gene activation of human coiled-coil domain containing 183                                                                                                                                                                                                                             |
| Symbol:       | CCDC183                                                                                                                                                                                                                                                                                                             |
| Components:   | GA113802G1, KIAA1984 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113802G2, KIAA1984 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113802G3, KIAA1984 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q5T5S1</a></u>                                                                                                                                                                                                                                                                                       |
| Synonyms:     | bA216L13.7; KIAA1984; PARF                                                                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)