

Product datasheet for **GA113821**

Human PIGY activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGY CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class Y
Symbol:	PIGY
UniProt ID:	Q3MUY2
Synonyms:	HPMRS6; PIG-Y
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)