

Product datasheet for **GA113831**

Human C18orf45 (TMEM241) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM241 CRISPRa kit - CRISPR gene activation of human transmembrane protein 241
Symbol:	TMEM241
UniProt ID:	Q24JQ0
Synonyms:	C18orf45; hVVT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)