


Product datasheet for **GA113847**

Human PPIL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPIL4 CRISPRa kit - CRISPR gene activation of human peptidylprolyl isomerase like 4
Symbol:	PPIL4
Components:	GA113847G1, PPIL4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113847G2, PPIL4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113847G3, PPIL4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WUA2
Synonyms:	HDCME13P
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)