

Product datasheet for **GA113848**

Human Progesterone Receptor Beta (PAQR8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAQR8 CRISPRa kit - CRISPR gene activation of human progesterone and adipoQ receptor family member 8
Symbol:	PAQR8
Components:	GA113848G1, Progesterone Receptor Beta gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113848G2, Progesterone Receptor Beta gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113848G3, Progesterone Receptor Beta gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8TEZ7
Synonyms:	C6orf33; LMPB1; MPRB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)