

Product datasheet for **GA113855**

Human SHISAL1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHISAL1 CRISPRa kit - CRISPR gene activation of human shisa like 1
Symbol:	SHISAL1
Components:	GA113855G1, SHISAL1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113855G2, SHISAL1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113855G3, SHISAL1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3SXP7</u>
Synonyms:	MGC125851; MGC125852
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)