

Product datasheet for **GA113864**

Human RIMBP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RIMBP3 CRISPRa kit – CRISPR gene activation of human RIMS binding protein 3
Symbol:	RIMBP3
Components:	GA113864G1, RIMBP3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113864G2, RIMBP3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113864G3, RIMBP3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UFD9</u>
Synonyms:	RIM-BP3; RIM-BP3.1; RIM-BP3.A; RIMBP3.1; RIMBP3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

