

Product datasheet for **GA113875**

Human Prostein (SLC45A3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC45A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 45 member 3
Symbol:	SLC45A3
UniProt ID:	Q96JT2
Synonyms:	IPCA-2; IPCA-6; IPCA-8; IPCA6; PCANAP2; PCANAP6; PCANAP8; PRST
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)