

Product datasheet for **GA113893**

Human TSPYL5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPYL5 CRISPRa kit - CRISPR gene activation of human TSPY like 5
Symbol:	TSPYL5
Components:	GA113893G1, TSPYL5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113893G2, TSPYL5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113893G3, TSPYL5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86VY4</u>
Synonyms:	KIAA1750
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)