

Product datasheet for **GA113911**

Human PSKH2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSKH2 CRISPRa kit - CRISPR gene activation of human protein serine kinase H2
Symbol:	PSKH2
Components:	GA113911G1, PSKH2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113911G2, PSKH2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113911G3, PSKH2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96QS6</u>
Synonyms:	protein serine kinase H2; serine/threonine kinase PSKH2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)