

Product datasheet for **GA113934**

Human PPP1R3F activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R3F CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 3F
Symbol:	PPP1R3F
Components:	GA113934G1, PPP1R3F gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113934G2, PPP1R3F gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113934G3, PPP1R3F gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZSY5</u>
Synonyms:	HB2E; LLOXNC01-7P3.1; R3F
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)