

Product datasheet for **GA113935**

Human **KCNK17** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | KCNK17 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 17 |
| Symbol: | KCNK17 |
| UniProt ID: | <u>Q96T54</u> |
| Synonyms: | K2p17.1; TALK-2; TALK2; TASK-4; TASK4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)