

Product datasheet for **GA113944**

Human SIGLECL1 (SIGLEC12) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SIGLEC12 CRISPRa kit - CRISPR gene activation of human sialic acid binding Ig like lectin 12 (gene/pseudogene)
Symbol:	SIGLEC12
UniProt ID:	Q96PQ1
Synonyms:	S2V; Siglec-XII; SIGLECL1; SLG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)