

Product datasheet for **GA113957**

Human **DYNC2I2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	WDR34 CRISPRa kit - CRISPR gene activation of human WD repeat domain 34
Symbol:	DYNC2I2
UniProt ID:	Q96EX3
Synonyms:	bA216B9.3; DIC5; FAP133; SRTD11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)