

## Product datasheet for **GA113960**

### Human BMERB1 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | BMERB1 CRISPRa kit - CRISPR gene activation of human bMERB domain containing 1  |
| Symbol:       | BMERB1  |
| Components:   | GA113960G1, BMERB1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113960G2, BMERB1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113960G3, BMERB1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q96MC5</a>  |
| Synonyms:     | FLJ32618  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)   |

#### Product images:



[View online »](#)