

Product datasheet for **GA113964**

Human KLC4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KLC4 CRISPRa kit - CRISPR gene activation of human kinesin light chain 4
Symbol:	KLC4
Components:	GA113964G1, KLC4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113964G2, KLC4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113964G3, KLC4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NSK0</u>
Synonyms:	bA387M24.3; KNSL8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)