

Product datasheet for **GA113965**

Human C9orf140 (SAPCD2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SAPCD2 CRISPRa kit - CRISPR gene activation of human suppressor APC domain containing 2
Symbol:	SAPCD2
Components:	GA113965G1, C9orf140 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113965G2, C9orf140 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113965G3, C9orf140 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UD0</u>
Synonyms:	C9orf140; p42.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)