


## Product datasheet for **GA113977**

### Human TMEM250 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM250 CRISPRa kit - CRISPR gene activation of human transmembrane protein 250
Symbol:	TMEM250
UniProt ID:	<u><a href="#">H0YL14</a></u>
Synonyms:	bA83N9.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)