

Product datasheet for **GA113978**

Human **TSR2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TSR2 CRISPRa kit – CRISPR gene activation of human TSR2 ribosome maturation factor
Symbol:	TSR2
Components:	GA113978G1, TSR2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113978G2, TSR2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113978G3, TSR2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q969E8
Synonyms:	DBA14; DTIP1A10; WGG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

