

## Product datasheet for **GA113979**

### Human **KCNH7** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | KCNH7 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily H member 7 |
| Symbol:       | KCNH7  |
| UniProt ID:   | <u><a href="#">Q9NS40</a></u>  |
| Synonyms:     | ERG3; HERG3; Kv11.3  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                  |

#### Product images:



[View online >](#)