

Product datasheet for **GA113999**

Human TGIF2LX activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TGIF2LX CRISPRa kit - CRISPR gene activation of human TGFB induced factor homeobox 2 like X-linked
Symbol:	TGIF2LX
Components:	GA113999G1, TGIF2LX gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113999G2, TGIF2LX gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113999G3, TGIF2LX gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IUE1</u>
Synonyms:	TGIFLX
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)