

Product datasheet for **GA114004**

Human EXOC3L2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EXOC3L2 CRISPRa kit - CRISPR gene activation of human exocyst complex component 3 like 2
Symbol:	EXOC3L2
Components:	GA114004G1, EXOC3L2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114004G2, EXOC3L2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114004G3, EXOC3L2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2M3D2</u>
Synonyms:	XTP7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)