

Product datasheet for **GA114009**

Human UNC119 binding protein C5orf30 (C5orf30) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C5orf30 CRISPRa kit - CRISPR gene activation of human chromosome 5 open reading frame 30
Symbol:	MACIR
Components:	GA114009G1, UNC119 binding protein C5orf30 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114009G2, UNC119 binding protein C5orf30 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114009G3, UNC119 binding protein C5orf30 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96GV9</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

