

## Product datasheet for **GA114015**

### Human TMEM41A activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM41A CRISPRa kit - CRISPR gene activation of human transmembrane protein 41A
Symbol:	TMEM41A
UniProt ID:	<a href="#">Q96HV5</a>
Synonyms:	2900010K02Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)