

Product datasheet for **GA114019**

Human TRAF4AF1 (KNSTRN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KNSTRN CRISPRa kit - CRISPR gene activation of human kinetochore localized astrin (SPAG5) binding protein
Symbol:	KNSTRN
UniProt ID:	Q9Y448
Synonyms:	C15orf23; HSD11; SKAP; TRAF4AF1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)