


## Product datasheet for **GA114024**

### Human **GADD45GIP1** activation kit by **CRISPRa**

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                           |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                              |
| Description:  | GADD45GIP1 CRISPRa kit - CRISPR gene activation of human GADD45G interacting protein 1                                                                                                                                                                                                                                    |
| Symbol:       | GADD45GIP1                                                                                                                                                                                                                                                                                                                |
| Components:   | GA114024G1, GADD45GIP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114024G2, GADD45GIP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114024G3, GADD45GIP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>Q8TAE8</u>                                                                                                                                                                                                                                                                                                             |
| Synonyms:     | CKBBP2; CKbetaBP2; CRIF1; MRP-L59; PLINP; PLINP-1; Plinp1; PRG6                                                                                                                                                                                                                                                           |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                   |

#### Product images:



[View online »](#)