

Product datasheet for **GA114030**

Human MLIP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MLIP CRISPRa kit - CRISPR gene activation of human muscular LMNA interacting protein
Symbol:	MLIP
Components:	GA114030G1, MLIP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114030G2, MLIP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114030G3, MLIP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5VWP3</u>
Synonyms:	C6orf142; CIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)