

Product datasheet for **GA114039**

Human ZNF700 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF700 CRISPRa kit - CRISPR gene activation of human zinc finger protein 700
Symbol:	ZNF700
Components:	GA114039G1, ZNF700 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114039G2, ZNF700 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114039G3, ZNF700 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H0M5</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

