

Product datasheet for **GA114040**

Human ZNF439 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF439 CRISPRa kit - CRISPR gene activation of human zinc finger protein 439
Symbol:	ZNF439
Components:	GA114040G1, ZNF439 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114040G2, ZNF439 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114040G3, ZNF439 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NDP4</u>
Synonyms:	DKFZp571K0837
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)