


Product datasheet for **GA114045**

Human ZFAND2A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZFAND2A CRISPRa kit - CRISPR gene activation of human zinc finger AN1-type containing 2A
Symbol:	ZFAND2A
Components:	GA114045G1, ZFAND2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114045G2, ZFAND2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114045G3, ZFAND2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8N6M9
Synonyms:	AIRAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)