


Product datasheet for **GA114048**

Human TGIF2LY activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TGIF2LY CRISPRa kit - CRISPR gene activation of human TGFB induced factor homeobox 2 like Y-linked
Symbol:	TGIF2LY
Components:	GA114048G1, TGIF2LY gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114048G2, TGIF2LY gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114048G3, TGIF2LY gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IUHQ</u>
Synonyms:	TGIFLY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)