

Product datasheet for **GA114054**

Human SEC11C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEC11C CRISPRa kit - CRISPR gene activation of human SEC11 homolog C, signal peptidase complex subunit
Symbol:	SEC11C
UniProt ID:	<u>Q9BY50</u>
Synonyms:	SEC11L3; SPC21; SPCS4C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)