

Product datasheet for **GA114068**

Human SPSB3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPSB3 CRISPRa kit - CRISPR gene activation of human splA/ryanodine receptor domain and SOCS box containing 3
Symbol:	SPSB3
UniProt ID:	Q6PJ21
Synonyms:	C16orf31; SSB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)