

## Product datasheet for **GA114070**

### Human C9orf123 (TMEM261) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	DMAC1 CRISPRa kit - CRISPR gene activation of human distal membrane arm assembly complex 1
Symbol:	DMAC1
UniProt ID:	<a href="#">Q96GE9</a>
Synonyms:	C9orf123
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)