

## Product datasheet for **GA114072**

### Human TRIM41 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | TRIM41 CRISPRa kit - CRISPR gene activation of human tripartite motif containing 41   |
| Symbol:       | TRIM41  |
| Components:   | GA114072G1, TRIM41 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114072G2, TRIM41 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA114072G3, TRIM41 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8WV44</a></u>   |
| Synonyms:     | RINCK   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)   |

#### Product images:



[View online »](#)