

Product datasheet for **GA114079**

Human FMNL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FMNL3 CRISPRa kit - CRISPR gene activation of human formin like 3
Symbol:	FMNL3
Components:	GA114079G1, FMNL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114079G2, FMNL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114079G3, FMNL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IVF7</u>
Synonyms:	FHOD3; FRL2; WBP-3; WBP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)