

Product datasheet for **GA114086**

Human TRIM47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRIM47 CRISPRa kit - CRISPR gene activation of human tripartite motif containing 47
Symbol:	TRIM47
Components:	GA114086G1, TRIM47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114086G2, TRIM47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114086G3, TRIM47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96LD4
Synonyms:	GOA; RNF100
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)