

Product datasheet for **GA114088**

Human L3MBTL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	L3MBTL4 CRISPRa kit - CRISPR gene activation of human L3MBTL histone methyl-lysine binding protein 4
Symbol:	L3MBTL4
UniProt ID:	Q8NA19
Synonyms:	HsT1031
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)