

Product datasheet for **GA114089**

Human **SLC25A46** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A46 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 46
Symbol:	SLC25A46
Components:	GA114089G1, SLC25A46 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114089G2, SLC25A46 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114089G3, SLC25A46 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96AG3</u>
Synonyms:	HMSN6B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)