

Product datasheet for **GA114101**

Human R3HDM4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	R3HDM4 CRISPRa kit - CRISPR gene activation of human R3H domain containing 4
Symbol:	R3HDM4
Components:	GA114101G1, R3HDM4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114101G2, R3HDM4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114101G3, R3HDM4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96D70</u>
Synonyms:	C19orf22
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)