

Product datasheet for **GA114106**

Human CDKN2AIPNL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CDKN2AIPNL CRISPRa kit - CRISPR gene activation of human CDKN2A interacting protein N-terminal like
Symbol:	CDKN2AIPNL
Components:	GA114106G1, CDKN2AIPNL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114106G2, CDKN2AIPNL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114106G3, CDKN2AIPNL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96HQ2</u>
Synonyms:	MGC13017
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)