

Product datasheet for **GA114113**

Human XRCC6BP1 (ATP23) activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: ATP23 CRISPRa kit - CRISPR gene activation of human ATP23 metalloproteinase and ATP synthase assembly factor homolog

Symbol: ATP23


Components: GA114113G1, XRCC6BP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)
GA114113G2, XRCC6BP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)
GA114113G3, XRCC6BP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)
1 CRISPRa-Enhancer vector, SKU GE100056
1 CRISPRa scramble vector, SKU GE100077

UniProt ID: [Q9Y6H3](#)

Synonyms: KUB3; XRCC6BP1

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)