


## Product datasheet for **GA114118**

### Human Protein kinase like protein SgK493 (PKDCC) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PKDCC CRISPRa kit - CRISPR gene activation of human protein kinase domain containing, cytoplasmic
Symbol:	PKDCC
UniProt ID:	<a href="#">Q504Y2</a>
Synonyms:	SGK493; Vlk
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)