

Product datasheet for **GA114127**

Human Plexin A4 (PLXNA4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLXNA4 CRISPRa kit - CRISPR gene activation of human plexin A4
Symbol:	PLXNA4
UniProt ID:	<u>Q9HCM2</u>
Synonyms:	FAYV2820; PLEXA4; PLXNA4A; PLXNA4B; PRO34003
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)