

Product datasheet for **GA114151**

Human ZNF804A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF804A CRISPRa kit - CRISPR gene activation of human zinc finger protein 804A
Symbol:	ZNF804A
Components:	GA114151G1, ZNF804A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA114151G2, ZNF804A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA114151G3, ZNF804A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7Z570
Synonyms:	C2orf10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)