

## Product datasheet for **GA114157**

### Human MYLK3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	MYLK3 CRISPRa kit - CRISPR gene activation of human myosin light chain kinase 3
Symbol:	MYLK3
UniProt ID:	<u><a href="#">Q32MK0</a></u>
Synonyms:	caMLCK; MLCK; MLCK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)